Cochlear™

Electrodo Slim Modiolar

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Referencias

1. Based on latest generation of cochlear implant electrodes manufactured by Cochlear available as at 1 July 2019.
2. Data on file - Investigation Report, 2017.
3. Aschendorff A, Briggs R, Brademann G, Helbig S, Hornung J, Lenarz T, Marx M, Ramos A, Stover T, Escudé B and James CJ. Clinical investigation of the Nucleus Slim Modiolar Electrode. *Audiology & Neurotology* 2017;22:169-179.
4. Shaul C, Dragovic AS, Stringer AK, O'Leary SJ, Briggs RJ. Scalar localisation of peri-modiolar electrodes and speech perception outcomes. *J Laryngol Otol*. 2018;132:1000-6.
5. Cuda D, Murri A. Cochlear implantation with the nucleus slim modiolar electrode (CI532): a preliminary experience. *Eur Arch Otorhinolaryngol*. 2017;274:4141-8.
6. Jiam N, Limb C. The impact of round window vs cochleostomy surgical approaches on interscalar excursions in the cochlea: Preliminary results from a flat-panel computed tomography study. *World Journal of Otorhinolaryngology-Head and Neck Surgery* (2016) 2, 142-147.
7. Ramos de Miguel A, Argudo AA, Borkoski Barreiro SA, Falcon Gonzalez JC, Ramos Macias A. Imaging evaluation of electrode placement and effect on electrode discrimination on different cochlear implant electrode arrays. *Eur Arch Otorhinolaryngol*. 2018 Jun;275(6):1385-1394.
8. Holder JT, Yawn RJ, Nassiri AM, Dwyer RT, Rivas A, Labadie RF, Gifford RH. Matched Cohort Comparison Indicates Superiority of Precurved Electrode Arrays. *Otol Neurotol* 2019 Oct;40(9):1160-1166.
9. Based on latest generation of cochlear implant electrodes manufactured by Cochlear, Advanced Bionics, MED-EL and Oticon Medical available as at 1 July 2019.
10. Briggs RJ, Tykocinski M, Laszig R, et al. Development and evaluation of the modiolar research array—multi-centre collaborative study in human temporal bones. *Cochlear Implants Int* 2011;12:129-39.
11. Holden LK, Finley CC, Firszt JB, Holden TA, Brenner C, Potts LG, et al. Factors affecting open-set word recognition in adults with cochlear implants. *Ear Hear* 2013; 34: 342-60.
12. Dowell R. Evidence about the effectiveness of cochlear implants for adults. Evidence based practice in audiology: evaluating interventions for children and adults with hearing impairment. Plural Publishing. 2012; 141-66.
13. Hey M, Wesarg T, Mewes A, Helbig S, Hornung J, Lenarz T, Briggs R, Marx M, Ramos A, Stöver T, Escudé B, James CJ, Aschendorff A. Objective, audiological and quality of life measures with the CI532 slim modiolar electrode. *Cochlear Implants Int*. 2019 Mar;20(2):80-90.
14. Friedmann DR, Kamen E, Choudhury B, Roland JT Jr. Surgical Experience and Early Outcomes With a Slim Perimodiolar Electrode. *Otol Neurotol*. 2019 40(3):e304-e310.
15. De Seta D, Torres R, Russo FY, Ferrary E, Kazmitcheff G, Heymann D, et al. Damage to inner ear structure during cochlear implantation: Correlation between insertion force and radio-histological findings in temporal bone specimens. *Hear Res* 2017; 344: 90-7.
16. Hughes ML, Abbas PJ. Electrophysiologic channel interaction, electrode pitch ranking, and behavioral threshold in straight versus perimodiolar cochlear implant electrode arrays. *J Acoust Soc Am* 2006; 119: 1538-47.
17. Xi X, Ji F, Han D, Hong M, Chen A. Electrode interaction in cochlear implant recipients: comparison of straight and contour electrode arrays. *ORL J Otorhinolaryngol Relat Spec* 2009; 71: 228-37.
18. Basta D, Todt I, Ernst A. Audiological outcome of the pull-back technique in cochlear implantation. *Laryngoscope* 2010; 120: 1391-6.
19. Saunders E, Cohen L, Aschendorff A, Shapiro W, Knight M, Stecker M, et al. Threshold, comfortable level and impedance changes as a function of electrode-modiolar distance. *Ear Hear* 2002; 23(1 Suppl): 285-405.
20. Parkinson AJ, Arcaroli J, Staller SJ, Arndt PL, Cosgriff A, Ebinger K. The nucleus 24 contour cochlear implant system: adult clinical trial results. *Ear Hear* 2002; 23(1 Suppl): 41S-8S.
21. Gomez Serrano M, Patel S, Harris R & Selvadurai D. Initial surgical and clinical experience with the Nucleus CI532 slim modiolar electrode in the UK. *Cochlear Implants International*, 2019; 20:4, 207-216.
22. Ramos-Macías A, Borkoski-Barreiro SA, Falcón-González JC, Ramos-de Miguel A. Hearing Preservation with the Slim Modiolar Electrode Nucleus CI532 Cochlear Implant: A Preliminary Experience. *Audiol Neurootol*. 2017;22:317-25.

Cochlear Ltd (ABN 96 002 618 073) 1 University Avenue, Macquarie University, NSW 2109, Australia Tel: +61 2 9428 6555 Fax: +61 2 9428 6352

Cochlear AG EMEA Headquarters, Peter Merian-Weg 4, 4052 Basel, Switzerland Tel: +41 61 205 8204 Fax: +41 61 205 8205

Cochlear Deutschland GmbH & Co. KG Karl-Wiechert-Allee 76A, 30625 Hannover, Germany Tel: +49 511 542 7750 Fax: +49 511 542 7770

Cochlear Europe Ltd 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey KT15 2HJ, United Kingdom Tel: +44 1932 26 3400 Fax: +44 1932 26 3426

Cochlear Austria GmbH CEE Office, Millennium Tower, 45th Floor, Handelskai 94-96, 1200 Vienna, Austria Tel: +43 1 37600 26 000

Cochlear Benelux NV Schaliënoevodreef 20 i, B-2800 Mechelen, Belgium Tel: +32 15 79 55 77

Cochlear Europe Limited Czech Branch Office, Kamenická 656/56, 170 00 Prague 7, Czech Republic Tel: +420 222 13 53 13

Cochlear France SAS 135 route de Saint Simon, CS 43574, 31035 Toulouse, France Tel: +33 5 34 63 85 85 (international) or 0805 200 016 (national) Fax: +33 5 34 63 85 80

Cochlear Italia S.r.l. Via Larga 33, 40138 Bologna, Italy Tel: +39 051 601 53 11 Fax: +39 051 39 20 62

Cochlear Middle East FZ-LLC Dubai Healthcare City, Al Razi Building 64, Block A, Ground Floor, Offices IR1 and IR2, Dubai, United Arab Emirates Tel: +971 4 818 4400 Fax: +971 4 361 8925

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© Cochlear Limited 2020. D1697906 V1 FEB20 Spanish Translation of D1668406 V3 FEB20
Cochlear Nordic AB Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden Tel: +46 31 335 14 61
Cochlear Nordic AB Menuetto Business Center, Huopalahdentie 24 A, 00350 Helsinki, Finland Tel: +358 20 735 0788